Gebrauchsanweisung Servomex

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Instructions for use Servomex

Artificial Nose (with Pressure Relief Valve)

REF 22353

1. Indication

The Servomex artificial nose is envisaged for spontaneously breathing patients, who breathe via an endotracheal tubus or a tracheal cannula. It can be used in clinical and in domestic applications.

2. Application instructions

The Servomex artificial nose is attached to a 15 mm standard connector of an endotracheal tubus, respectively a tracheal cannula. Here it is important to ensure that it is firmly attached!

The Servomex artificial nose is designed for a breathing volume of 50 – 1000 ml.

If an additional supply of oxygen is required, connect the oxygen connector via an oxygen tube to an oxygen source (oxygen concentrator, oxygen wall connection, etc.).

The normal moistening effect begins after a few breaths.

The Servomex artificial nose is not fitted with any metal parts and does not conduct electricity.

3. Contraindications

Do not use the Servomex artificial nose with other than the recommended breathing volume. If the breathing volume is too low, the dead air space increases, which can lead to an increase in CO₂. If the breathing volume is too great, this leads to insufficient moistening of the air.

We also recommend that the device should not be used for the following patients:

- Patients with severe exsiccosis
- Patients with excessive secretion in the respiratory tracts
- Patients with severe air leak current (the air that is breathed out does not flow through the artificial nose)
- Patients who use a speech valve

4. Precautionary measures

The Servomex artificial nose should be replaced at the latest after 24 hours. Depending on requirements, it can be changed more frequently.

If there is an increase in the total pulmonary resistance, for instance through increased secretory accumulation in the artificial nose, the artificial nose should be removed.

The Servomex artificial nose is designed for being used only once!

Please immediately dispose of the artificial nose when it has been replaced!

5. Technical data (DIN EN ISO 9360 - 2:2002)

Recommended breathing volumes:	50 – 1000 ml
Internal volume:	17 ml
Moisture loss:	VT = 250 ml : 8,5 mg/ VT = 500 ml : 10.6 mg/l VT = 750 ml : 13.0 mg/l VT = 1000ml : 17.0 mg/l
Drop in pressure: Inspiration:	0.5 l / s = 0,08 h Pa 1.0 l / s = 0,13 h Pa 1.5 l / s = 0,16 h Pa
Expiration:	0.5 l / s = 0,09 h Pa 1.0 l / s = 0,14 h Pa 1.5 l / s = 0,17 h Pa
Weight:	6 g
Length:	29 mm
Connections:	– 15 mm – for attaching to a 15 mm standard connector of an endotracheal tubus, respectively a tracheal cannula – 5 mm plug junction for oxygen inlet – Safety valve for forced breathing in/breathing out Secretion suction
Safety valve:	Enables high short-term tidal volumes and facilitates coughing out great volumes of secrete
Environmental compatibility:	When the product is incinerated this produces hydrocarbon, carbon dioxide, water and carbon oxide.

6. Sign and explanation for:

Disposable product

CE – symbol und number of the named place

Follow Operating instructions

Batch number/ Date of manufacture

Usable until

servona
Uns verbindet Gesundheit.

Servona GmbH

Biberweg 24-26, D-53842 Troisdorf
Tel. +49(0)22 41.93 22-0
Fax +49(0)22 41.93 22-277
info@servona.com | www.servona.com